

The logo for LINC (Liaison for Interventional Neurovascular Care) features a stylized, colorful graphic of a brain or vessel with red and orange highlights, set against a dark blue background. The letters "LINC" are positioned to the right of this graphic.

LINC

12-month Results from the AcoArtII BTK Study

Xin Jia Wei Guo

Vascular and endovascular surgery center
Chinese PLA General Hospital, Beijing

on behalf of the AcoArt II Investigators

Disclosure


Speaker name: Xin Jia

I have the following potential conflicts of interest to report:

- Consulting
 - Employment in industry
 - Stockholder of a healthcare company
 - Owner of a healthcare company
 - Other(s)
-
- I do not have any potential conflict of interest

AcoArt II Study Overview

- Primary Endpoint
 - 6-mon Angiographic Primary patency (PP) – freedom from target lesion occlusion, CD-TLR and major amputation.
- Key Secondary Endpoint
 - Late lumen loss (LLL) at 6mon
 - CD-TLR at 6mon, 12mon
 - Ulcer healing Rate at 12mon
 - Major amputation at 12mon
- Safety Endpoint
 - MAE (death, major amputation, CD-TLR) at 30 days, 6mon, 12mon



Study Device - Litos and Tulip DCB

- ***Litos (0.014") & Tulip (0.018")***® produced by Acotec Scientific Co., Ltd, Beijing, China
- Paclitaxel dose:
 - 3.3µg/mm²
- Excipient: Mg-stearate
 - Durable paclitaxel carrier for optimized drug delivery
 - Uniform coating of balloon treatment area
- Balloon sizes
 - Diameter(mm):
 - *Litos*: 2.0-4.0, with 0.5mm increment
 - *Tulip*: 2.0 -4.0, with 0.5mm increment
 - Length(mm):
 - 20, 30, 40, 60, 80, 100, 120, 150, 200, 250, 300, treat long lesion easily



AcoArt II Baseline

	DCB (N=61)	PTA (N=59)	P Value
Age, yrs	70.7±7.4	70.8±9.0	0.95
Male (%)	59%(36)	61%(36)	0.82
History of Risk Factors			
CHD	36%(22)	34%(20)	0.81
Hypertension	82%(50)	75%(44)	0.33
Hyperlipidemia	41%(25)	27%(16)	0.11
Diabetes	74%(45)	71%(42)	0.75
Current smoker	26%(16)	27%(16)	0.60
Current alcoholic	16%(10)	10%(6)	0.32
Rutherford			
3	2%(1)	0	
4	44%(27)	41%(24)	0.76
5	39%(24)	47%(28)	
6	15%(9)	12%(7)	
ABI*	0.56±0.27	0.51±0.31	0.41

Values are % (n/N), or mean ±SD. *Calculated by excluding ABI≥1.3.

Target Lesion Characteristics



	DCB	PTA	P
Patients (No)	61	59	
Lesions (No)	65	66	
Target vessels			
TA	1.5%(1)	1.5%(1)	
TA+PTA	9.2%(6)	6.1%(4)	
TA+PA	4.6%(3)	7.6%(5)	0.59
ATA	49%(32)	50%(33)	
PTA	26%(17)	18%(12)	
PA	9.2%(6)	17%(11)	
Target lesion length (mm)	177±86	186±82	0.56
Minimal lumen diameter (mm)	0.20±0.41	0.15±0.33	0.43
Reference vessel diameter (mm)	2.55±0.34	2.53±0.35	0.71
Diameter stenosis(%)	92%	94%	0.43
CTO (%)	75.4% (49/65)	78.8% (52/66)	0.54

Values are % (n/N), or mean ±SD.

TA= tibiofibular artery; ATA= anterior tibial artery; PTA= posterior tibial artery; PA= peroneal artery

Procedural Data

	DCB	PTA	P
Treatment balloon diameter (mm)	2.72±0.25	2.64±0.25	0.07
Provisional stenting %	0	1.5% (1/66)	0.32
Device success	100%	100%	/




Appropriate diameter of DCB

- 3.0-3.5mm for upper calf
- 3.0-2.5mm for middle calf
- 2.5mm for ankle area


6-month results - Primary Endpoint

	DCB (N=48)	PTA (N=46)	P
Primary Patency	77.1% (37/48)	28.3% (13/46)	
Occlusion	8	27	<0.001
TLR	3	12	
Amputation	1	1	

• Primary patency was defined as *Freedom from occlusion and CD-TLR and major amputation above ankle .*




Primary Endpoint (KM 6 month)




12-month Results

	DCB	PTA	<i>P-value</i>
<i>CD-TLR</i>	8.5%(5/59)	23.2%(13/56)	0.030
<i>Ulcer healing rate</i>	83.9%(26/31)	75.0%(24/32)	0.393
<i>Major Amputation</i>	1.7%(1/59)	1.8%(1/56)	0.971



KM 12-month Freedom from CD-TLR




Major adverse event (MAE) - 12mon

	DCB	PTA	P value
Death	1.7% (1/59)	3.6% (2/56)	0.532
Major amputation	1.7% (1/59)	1.8% (1/56)	0.971
CD-TLR	8.5% (5/59)	23.2% (13/56)	0.030

Primary Endpoint Comparison

	AcoArt II	Lutonix BTK
Age, yrs	70.7	72.9
RVD(mm)	2.53	2.54
Lesion Length(cm)	18.2	10.5
DS(%)	93%	86.0%
CTO(%)	77%	35%
Paclitaxel dose	3.3ug/mm ²	2.0ug/mm ²

Primary Patency@6 mon



Compared with Lutonix BTK, AcoArt II

- More complex lesions with longer length and higher occlusion rate
- More significant difference between DCB and PTA arm in primary patency
- Use the DCB with higher paclitaxel dose

Conclusion

- AcoArt II demonstrated the safety and efficacy of Tulip/Litos DCB in treating BTK lesions.
- Further add to the evidence of DCB for BTK, which encourages utilization of DCB for future.
- More well designed study based on BTK-dedicated DCB may provide further data.

The logo for LINC (Lancet Interdisciplinary Network for Cardiovascular Care) features a stylized, colorful graphic of a heart or vessel in shades of red, orange, and yellow, set against a dark blue background. The letters "LINC" are positioned to the right of this graphic.

LINC

12-month Results from the AcoArtII BTK Study

Xin Jia Wei Guo

Vascular and endovascular surgery center
Chinese PLA General Hospital, Beijing

on behalf of the AcoArt II Investigators