

***The Magic Touch - PTA Sirolimus
Coated Balloon - Insights from
pre-clinical evaluation***

Aloke Finn

CVPath Institute, Inc.

Gaithersburg, MD, USA

Disclosure

Within the past 12 months, I or my spouse/partner have had a financial interest/arrangement or affiliation with the organization(s) listed below.

Employment in industry: No

Honorarium:

Amgen; Abbott Vascular; Biosensors; Boston Scientific; Celonova; Cook Medical; CSI; Lutonix Bard; Sinomed; Terumo Corporation.

Institutional grant/research support:

R01 HL141425 Leducq Foundation Grant; 480 Biomedical; 4C Medical; 4Tech; Abbott; Accumedical; Amgen; Biosensors; Boston Scientific; Cardiac Implants; Celonova; Claret; Concept Medical; Cook; CSI; DuNing; Edwards; Emboline; Endotronix; Envision Scientific; Lutonix/Bard; Gateway; Lifetech; Limflo; MedAlliance; Medtronic; Mercator; Merrill; Microport; Microvention; Mitraalign; Mitraassist; NAMSA; Nanova; Neovasc; NIPRO; Novogate; Occulotech; Orbus Neich; Phenox; Profusa; Protembis; Qool; Recor; Senseonics; Shockwave; Sinomed; Spectranetics; Surmodics; Symic; Vesper; W.L. Gore; Xeltis.

Owner of a healthcare company: No

Stockholder of a healthcare company: No


Elements of an Effective DCB Formulation

- ✓ Must deliver large quantities of the drug within seconds
- ✓ Distribute within the media in the first few days
- ✓ Therapeutic drug levels must be maintained for more than 4 weeks
- ✓ Must allow rapid healing as compared to DES
- ✓ No need for long-term anti-platelet therapy Biologic effects must be observed by histology at 28-days
- ✓ Effective drug delivery to target tissue while avoiding non-target effect

Drug Coated Balloon Devices (Peripheral artery)


Device	Company	Coating	Drug dose ($\mu\text{g}/\text{mm}^2$)	CE mark*
Advance 18 PTX™	Cook Medical, Bloomington, IN, USA	Paclitaxel	3.0	Yes
Cotavance®	Bayer Schering Pharma AG, Berlin, Germany	Paclitaxel–iopromide	3.0	Yes
Freeway™	Eurocor, Bonn, Germany	Paclitaxel–shellac	3.0	Yes
In.Pact™ Admiral,	Medtronic Vascular, Santa Clara, CA, USA	Paclitaxel–urea	3.5	Yes
Lutonix® 035 DCB	BARD, Murray Hill, NJ, USA	Paclitaxel–polysorbate/sorbitol	2.0	Yes
Ranger	Boston Scientific	Paclitaxel–Acetyl Tributyl Citrate	2.0	Yes
Passeo-18 Lux®	Biotronik, Bülach, Switzerland	Paclitaxel–butyryl-tri-hexyl citrate	3.0	No → Yes
Stellarex®	Covidien, Mansfield, MA, USA	Paclitaxel	2.0	Yes
SurVeil™DCB	SurModics, MN, USA	Paclitaxel-proprietary photolink®	2.0	No → No

Byrne RA, Joner M. et al. Nat Rev Cardiol. 2014;11:13-23




Sirolimus and Paclitaxel

Sirolimus




Solid line : FKBP12 binding site
Dotted line : mTOR binding site

Sirolimus inhibits mTOR and is part of the phosphatidylinositol kinase-related family of serine/threonine kinases.



Habib et al. *Interventional Cardiology Clinics* 2016;5:321-329


Paclitaxel



Paclitaxel binds to beta-tubulin and impairs microtubular disassembly and halts the cell cycle between G₂ and M

Gupta ML et al. *PNAS* 2003;100:6394-6397

Dose dependent effect and toxicity of paclitaxel eluting stent in porcine coronary artery (28days)




BMS

**Paclitaxel
0.2µg/stent**


**Paclitaxel
15µg/stent**

**Paclitaxel
187µg/stent**


Neointimal area



(%) % neointimal fibrin



**(%) % medial wall cell
necrosis**




Sirolimus offers potential benefits over Paclitaxel

Common anti-proliferative drug for DCB is currently **PACLITAXEL**, however, **SIROLIMUS** (rapamycin) offers potential benefits over Paclitaxel.

	SIROLIMUS (OR ANALOGS)	PACLITAXEL
Inhibition of SMC proliferation	++	++
Inhibition of SMC migration	++	+
Inhibition of EC proliferation	++	++
Pro-apoptotic effects	(+)	++
Therapeutic range	WIDE	NARROW
Safety margin	10'000 fold	100 fold
Anti-Restenotic impact	++	+
Anti-inflammatory properties	++	(+) / -
Tissue Absorption	SLOW	FAST
Tissue Retention	SHORT	LONG

Risk of Death following Application of PES and PCB in Femoropopliteal artery

Random effects forest plot of all-cause death at 2 years




Meta-Analysis of RCT

Causes of Death

	Paclitaxel-Coated Balloon (IN.PACT SFA) at 3 Years ^{10,62}		Paclitaxel-Coated Stent (ZILVER PTX) at 2 Years ^{19,23}	
	Paclitaxel	Control	Paclitaxel	Control
Cardiovascular	9	0	18	8
Cancer	2	2		
Infectious	5	0		
Pulmonary	3	0		
Other	3	0	NA	NA

Random effects forest plot of all-cause death at 4 to 5 years



Katsanos K, et al. J Am Heart Assoc. 2018;7:e011245

Sirolimus Coated Balloons – Technical challenges

- **Enhance tissue absorption**

- Difficult to get sirolimus to enter into arterial tissue within 30 to 180 seconds of balloon dilatation; hence some kind of “instant glue” is required to transfer the drug from the balloon to the tissue efficiently


- **Extend tissue retention**

- Sirolimus must be continuously delivered over time, so some form of “time release mechanism” must be employed to maintain therapeutic levels

MAGIC TOUCH – Sirolimus Coated Balloon

- MAGICTOUCH® – SCB is Sirolimus Coated Balloon to treat coronary artery disease
- Delivers drug in 60 seconds
- Sub-micron particles
- Supersedes limitations of Paclitaxel
- No permanent metallic cage
- No more toxic effects
- Address many limitations of DES

Nothing Leaves Behind



Pre-clinical study 1; swine iliofemoral artery

Day 0


DCB treatment in
iliofemoral arteries

Day 28

Euthanize



Yorkshire domestic
Swine n=2
(total 4 vessels)




✓ Collect iliofemoral
arteries and
histopathologic
analysis

Control balloon

XtremeTouch SCB 3x





Pre-procedure



Follow-up 28day





Pre-procedure




Follow-up 28day

Histopathologic analysis for SCB treated site


Control balloon




XtremeTouch SCB 3x




Medial proteoglycan collagen score



Medial SMC loss score (Depth)



Medial SMC loss score (Circumf)



Pre-clinical study 2; swine coronary ISR lesion

Day 0

Vessel injury and
BMS implantation
in coronary arteries

Yorkshire
domestic
Swine n=7
(total 16 vessels)



Day 28

Angiogram and
treatment of ISR with
DCB or plain balloon

Magic Touch 1x,
Magic Touch 3x,
POBA 1x,
POBA 3x

Day 56

Euthanize

Collect

- ✓ Coronary artery
- ✓ Blood
- ✓ Tissue and Organs



Pharmacokinetics
and
Histopathology
Analysis

Pharmacokinetic data

Blood concentration [Sirolimus, ng/mL] in DCB 3x animal

Baseline	15min	1h	24h	33day
BQL	23.3±0.35	17.0±1.20	1.11±0.15	BQL

BQL: Below Quantitation Limit (0.250 ng/ml)

Artery and surrounding myocardium [Sirolimus, ng/g] in DCB treated animal (33-day after DCB)

Animal #	Treatment	Artery			Surrounding myocardium	
		DCB treated section	Proximal section	Distal section	Middle of treated area	Distal of treated area
Animal #1	DCB 3x (LAD)	525	BQL	BQL	BQL	BQL
Animal #2	DCB1x (RCA)	118	BQL	BQL	BQL	BQL

BQL: 50.0 ng/g

Myocardium and other organs [Sirolimus, ng/g] in DCB treated animal (33-day after DCB)

Animal #	Treatment	AA	AL	AP	AS	APRV	MA	ML	MP	MS	MPRV
Animal #1	DCB 3x (LAD)	BQL	BQL	BQL	BQL	BQL	BQL	BQL	BQL	BQL	BQL
Animal #2	DCB1x (RCA)	BQL	BQL	BQL	BQL	BQL	BQL	BQL	BQL	BQL	BQL

Animal #	Treatment	Lung	Liver	Bone Marrow	Spleen	Kidney
Animal #1	DCB 3x (LAD)	BQL	2.76	BQL	BQL	BQL
Animal #2	DCB1x (RCA)	BQL	BQL	BQL	BQL	BQL

AA: apical anterior, AS: apical septum, AP: apical posterior, AL: apical lateral, ARV: apical right ventricle
MA: mid anterior, MS: mid septum, MP: mid posterior, ML: mid lateral, MRV: mid right ventricle


BQL: 2.5 ng/g

Histopathologic analysis for SCB vs uncoated POBA treated site


MagicTouch SCB 1x

MagicTouch SCB 3x


Control POBA 3x




Medial PGs/Collagen score



Medial SMC loss score



Neointimal fibrin score







* All section base analysis

New neointimal growth was suppressed in SCB treatment; Major IEL disruption lesion

MagicTouch SCB 3x

Control POBA 3x



* All section base analysis

Summary

- ✓ Neointimal hyperplasia is the major cause of restenosis. Pharmacologic agents used on DES and DCB mainly inhibit SMC proliferation.
- ✓ The most commonly used anti-proliferative agent for DCB is currently Paclitaxel (Ptx). However cytotoxic Ptx has narrow therapeutic window compared with sirolimus and its analogues.
- ✓ Recent meta-analysis suggested that DES and DCB using PTX associated with higher death rates at 2 to 5 years, although the mechanisms and results still remain controversial.
- ✓ MagicTouch Sirolimus DCB has novel coating technology allowing better tissue retention of sirolimus drug.
- ✓ Our pre-clinical studies support efficacy of MagicTouch Sirolimus Coated Balloons.
- ✓ Anti-restenotic effect with minimal non-target organ drug levels were observed.
- ✓ MagicTouch DCB will be an important non-paclitaxel option for peripheral and coronary intervention

Acknowledgments

CVPath Institute

Atsushi Sakamoto, MD, PhD

Masayuki Mori, MD, PhD

Yu Sato, MD

Rika Kawakami, MD, PhD

Hiroyuki Jinnouchi, MD

Anne Cornelissen, MD

Liang Guo, PhD

Robert Kutyz, MS

Russ Jones


Abebe Atiso, HT

Jinky Beyer

Lila Adams, HT

Frank D Kolodgie, PhD


Renu Virmani, MD



≥30 days Outcomes after Paclitaxel-DES Implantation in Coronary Artery: CVPath autopsy registry

Late stent thrombosis; 69-year old man (SVG)

Late stent thrombosis; 48-year old man (LAD)



Patients, n	66
Age, yrs	60 ± 12
Male	61 (77%)
Stent related cardiac death	22 (33%)

Lesions, n	94
Duration, days	211 (118-383)
Stent length, mm	21.5 (16.0-30.0)
Stent outcome: Thrombosis	<u>25 (27%)</u>
Restenosis	5 (5%)
<u>Struts with fibrin, %</u>	<u>50 (36-69)</u>
<u>Uncovered struts, %</u>	<u>20 (8-49)</u>
<u>Malapposition, n(%)</u>	<u>18 (19%)</u>

Cause of late thrombosis in PES

Late thrombosis (≥ 30days), n	25
AMI, penetration/prolapse of necrotic core	5 (20%)
Bifurcation	9 (36%)
Long/overlapping stents	2 (0.1%)
Underexpansion	1 (0.04%)
Isolated uncovered struts	1 (0.04%)
Localized hypersensitivity reaction	0 (0.0%)
Malapposition from excessive to fibrin *	7 (28%)

*We have no late stent thrombosis case of sirolimus-eluting stent due to excessive fibrin around stent strut.

***The Magic Touch - PTA Sirolimus
Coated Balloon - Insights from
pre-clinical evaluation***

Aloke Finn

CVPath Institute, Inc.

Gaithersburg, MD, USA